

Interleukin-34

Catalog No: #AP84173



Package Size: #AP84173-1 50ug #AP84173-2 100ug #AP84173-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-34
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-34,C16orf77
Accession No.	Uniprot ID:Q6ZMJ4
Uniprot	Q6ZMJ4
GeneID	146433;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL34

Note: This product is for in vitro research use only