

Interleukin-6 receptor subunit alpha

Catalog No: #AP84174



Package Size: #AP84174-1 50ug #AP84174-2 100ug #AP84174-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-6 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IL-6 receptor subunit alpha,IL-6R 1,CD126,Il6ra
Accession No.	Uniprot ID:P22273
Uniprot	P22273
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Il6r

Note: This product is for in vitro research use only