

Golgi membrane protein 1

Catalog No: #AP84181



Package Size: #AP84181-1 50ug #AP84181-2 100ug #AP84181-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Golgi membrane protein 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Golgi membrane protein GP73,Golgi phosphoprotein 2,Golph2 |
| Accession No. | Uniprot ID:Q91XA2 |
| Uniprot | Q91XA2 |
| GeneID | 105348; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Golm1

Note: This product is for in vitro research use only