

## NF-kappa-B inhibitor beta

Catalog No: #AP84187



Package Size: #AP84187-1 50ug #AP84187-2 100ug #AP84187-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NF-kappa-B inhibitor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NF-kappa-BIB,I-kappa-B-beta,IkB-B,Ikbb
Accession No.	Uniprot ID:Q60778
Uniprot	Q60778
GeneID	18036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Nfkbib

Note: This product is for in vitro research use only