

Hemoglobin subunit alpha-1/2

Catalog No: #AP84191



Package Size: #AP84191-1 50ug #AP84191-2 100ug #AP84191-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hemoglobin subunit alpha-1/2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-1/2-globin,Hemoglobin alpha-1/2 chain,Hba-a1
Accession No.	Uniprot ID:P01946
Uniprot	P01946
GeneID	25632;360504;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hba1

Note: This product is for in vitro research use only