

Glycophorin-A

Catalog No: #AP84200



Package Size: #AP84200-1 50ug #AP84200-2 100ug #AP84200-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glycophorin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MN sialoglycoprotein,PAS-2,Sialoglycoprotein alpha,CD235a,GPA
Accession No.	Uniprot ID:P02724
Uniprot	P02724
GeneID	2993;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GYPA

Note: This product is for in vitro research use only