

Heat shock 70 kDa protein 4

Catalog No: #AP84201



Package Size: #AP84201-1 50ug #AP84201-2 100ug #AP84201-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock 70 kDa protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ischemia responsive 94 kDa protein,Irp94
Accession No.	Uniprot ID:O88600
Uniprot	O88600
GeneID	266759;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspa4

Note: This product is for in vitro research use only