

Heat shock protein beta-1

Catalog No: #AP84202



Package Size: #AP84202-1 50ug #AP84202-2 100ug #AP84202-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock protein beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	HspB1,Heat shock 27 kDa protein,HSP 27,HSP27
Accession No.	Uniprot ID:Q5S1U1
Uniprot	Q5S1U1
GeneID	493184;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSPB1

Note: This product is for in vitro research use only