

Insulin-like growth factor-binding protein 5

Catalog No: #AP84204



Package Size: #AP84204-1 50ug #AP84204-2 100ug #AP84204-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor-binding protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IBP-5,Igfbp-5
Accession No.	Uniprot ID:P24594
Uniprot	P24594
GeneID	25285;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Igfbp5

Note: This product is for in vitro research use only