

## Hepatitis A virus cellular receptor 2 homolog

Catalog No: #AP84205



Package Size: #AP84205-1 50ug #AP84205-2 100ug #AP84205-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hepatitis A virus cellular receptor 2 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HAVcr-2,T-cell immunoglobulin and mucin domain-containing protein 3,TIMD-3,T-cell immunoglobulin mucin receptor 3,TIM-3,T-cell membrane protein 3,Tim3,Timd3
Accession No.	Uniprot ID:Q8VIM0
Uniprot	Q8VIM0
GeneID	171285;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Havcr2

Note: This product is for in vitro research use only