

## Inhibitor of nuclear factor kappa-B kinase subunit beta

Catalog No: #AP84207



Package Size: #AP84207-1 50ug #AP84207-2 100ug #AP84207-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Inhibitor of nuclear factor kappa-B kinase subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	I-kappa-B-kinase beta,I-kappa-B kinase 2,IKK2,Nuclear factor NF-kappa-B inhibitor kinase beta,NFKBKB
Accession No.	Uniprot ID:Q95KV0
Uniprot	Q95KV0
GeneID	281854;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IKBKB

Note: This product is for in vitro research use only