

Interleukin-4 receptor subunit alpha

Catalog No: #AP84209



Package Size: #AP84209-1 50ug #AP84209-2 100ug #AP84209-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-4 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IL-4 receptor subunit alpha,CD124,IL4ra
Accession No.	Uniprot ID:Q63257
Uniprot	Q63257
GeneID	25084;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL4r

Note: This product is for in vitro research use only