

Kinesin heavy chain isoform 5A

Catalog No: #AP84211



Package Size: #AP84211-1 50ug #AP84211-2 100ug #AP84211-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kinesin heavy chain isoform 5A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Kinesin heavy chain neuron-specific 1, Neuronal kinesin heavy chain, NKHC, Kiaa4086, Kif5, Nkhc1
Accession No.	Uniprot ID:P33175
Uniprot	P33175
GeneID	16572;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Kif5a

Note: This product is for in vitro research use only