

Hydroxymethylglutaryl-CoA synthase, cytoplasmic

Catalog No: #AP84226



Package Size: #AP84226-1 50ug #AP84226-2 100ug #AP84226-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hydroxymethylglutaryl-CoA synthase, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HMG-CoA synthase,3-hydroxy-3-methylglutaryl coenzyme A synthase,HMGCS
Accession No.	Uniprot ID:Q01581
Uniprot	Q01581
GeneID	3157;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HMGCS1

Note: This product is for in vitro research use only