

## Interleukin-3

Catalog No: #AP84232



Package Size: #AP84232-1 50ug #AP84232-2 100ug #AP84232-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IL-3,Hematopoietic growth factor,Mast cell growth factor,MCGF,Multipotential colony-stimulating factor,P-cell-stimulating factor,IL-3
Accession No.	Uniprot ID:P04823
Uniprot	P04823
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL3

Note: This product is for in vitro research use only