

Interleukin-6

Catalog No: #AP84233



Package Size: #AP84233-1 50ug #AP84233-2 100ug #AP84233-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	ChIL-6
Accession No.	Uniprot ID:Q90Y10
Uniprot	Q90Y10
GeneID	395337;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL6

Note: This product is for in vitro research use only