

Intercellular adhesion molecule 3

Catalog No: #AP84236



Package Size: #AP84236-1 50ug #AP84236-2 100ug #AP84236-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Intercellular adhesion molecule 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ICAM-3,CD50
Accession No.	Uniprot ID:Q28125
Uniprot	Q28125
GeneID	281840;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ICAM3

Note: This product is for in vitro research use only