

Hepatoma-derived growth factor

Catalog No: #AP84237



Package Size: #AP84237-1 50ug #AP84237-2 100ug #AP84237-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------|
| Product Name | Hepatoma-derived growth factor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | HDGF,Tdrm1 |
| Accession No. | Uniprot ID:P51859 |
| Uniprot | P51859 |
| GeneID | 15191; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hdgf

Note: This product is for in vitro research use only