

Hepatocyte growth factor-like protein

Catalog No: #AP84238



Package Size: #AP84238-1 50ug #AP84238-2 100ug #AP84238-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hepatocyte growth factor-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Macrophage stimulatory protein,Macrophage-stimulating protein,MSP,HGFL
Accession No.	Uniprot ID:Q24K22
Uniprot	Q24K22
GeneID	514667;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MST1

Note: This product is for in vitro research use only