

High mobility group protein HMG-I/HMG-Y

Catalog No: #AP84239



Package Size: #AP84239-1 50ug #AP84239-2 100ug #AP84239-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	High mobility group protein HMG-I/HMG-Y
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	HMG-I(Y),High mobility group AT-hook protein 1,High mobility group protein A1,HMGIY
Accession No.	Uniprot ID:Q6URC2
Uniprot	Q6URC2
GeneID	442946;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HMGA1

Note: This product is for in vitro research use only