

Hematopoietic prostaglandin D synthase

Catalog No: #AP84240



Package Size: #AP84240-1 50ug #AP84240-2 100ug #AP84240-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Hematopoietic prostaglandin D synthase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | H-PGDS,GST class-sigma,Glutathione S-transferase,Glutathione-dependent PGD synthase,Glutathione-requiring prostaglandin D synthase,Prostaglandin-H2 D-isomerase,Gsts,Pgds,Ptgds2 |
| Accession No. | Uniprot ID:Q9JHF7 |
| Uniprot | Q9JHF7 |
| GeneID | 54486; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hpgds

Note: This product is for in vitro research use only