

Heat shock protein HSP 90-alpha

Catalog No: #AP84241



Package Size: #AP84241-1 50ug #AP84241-2 100ug #AP84241-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock protein HSP 90-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Heat shock 86 kDa,HSP 86,Tumor-specific transplantation 86 kDa antigen,TSTA,Hsp86,Hsp86-1,Hspca
Accession No.	Uniprot ID:P07901
Uniprot	P07901
GeneID	15519;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hsp90aa1

Note: This product is for in vitro research use only