

Inter-alpha-trypsin inhibitor heavy chain H4

Catalog No: #AP84242

Package Size: #AP84242-1 50ug #AP84242-2 100ug #AP84242-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inter-alpha-trypsin inhibitor heavy chain H4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ITI heavy chain H4,Inter-alpha-trypsin inhibitor family heavy chain-related protein,IHRP,Plasma kallikrein sensitive glycoprotein 120,Gp120,IHRP,ITIHL1,PK120,PRO1851
Accession No.	Uniprot ID:Q14624
Uniprot	Q14624
GeneID	3700;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITIHL4

Note: This product is for in vitro research use only