

Tyrosine-protein kinase JAK2

Catalog No: #AP84243



Package Size: #AP84243-1 50ug #AP84243-2 100ug #AP84243-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tyrosine-protein kinase JAK2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Janus kinase 2,JAK-2
Accession No.	Uniprot ID:Q62120
Uniprot	Q62120
GeneID	16452;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Jak2

Note: This product is for in vitro research use only