

Tyrosine-protein kinase JAK2

Catalog No: #AP84244



Package Size: #AP84244-1 50ug #AP84244-2 100ug #AP84244-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tyrosine-protein kinase JAK2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Janus kinase 2,JAK-2
Accession No.	Uniprot ID:O19064
Uniprot	O19064
GeneID	397201;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:JAK2

Note: This product is for in vitro research use only