

Interleukin-12 subunit beta

Catalog No: #AP84246



Package Size: #AP84246-1 50ug #AP84246-2 100ug #AP84246-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-12 subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	IL-12B,Cytotoxic lymphocyte maturation factor 40 kDa subunit,CLMF p40,IL-12 subunit p40
Accession No.	Uniprot ID:Q28268
Uniprot	Q28268
GeneID	403976;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL12B

Note: This product is for in vitro research use only