

Interleukin-2 receptor subunit beta

Catalog No: #AP84248



Package Size: #AP84248-1 50ug #AP84248-2 100ug #AP84248-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Interleukin-2 receptor subunit beta |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | IL-2 receptor subunit beta,High affinity IL-2 receptor subunit beta,p70-75,CD122 |
| Accession No. | Uniprot ID:P26896 |
| Uniprot | P26896 |
| GeneID | 25746; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL2rb

Note: This product is for in vitro research use only