

Interleukin-7

Catalog No: #AP84249



Package Size: #AP84249-1 50ug #AP84249-2 100ug #AP84249-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	IL-7
Accession No.	Uniprot ID:Q9N2G6
Uniprot	Q9N2G6
GeneID	397253;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL7

Note: This product is for in vitro research use only