

Inosine triphosphate pyrophosphatase

Catalog No: #AP84251



Package Size: #AP84251-1 50ug #AP84251-2 100ug #AP84251-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inosine triphosphate pyrophosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ITPase, Non-canonical purine NTP pyrophosphatase, Non-standard purine NTP pyrophosphatase, Nucleoside-triphosphate diphosphatase, Nucleoside-triphosphate pyrophosphatase, NTPase, Putative oncogene protein hlc14-06-p, C20orf37, My049, OK/SW-cl.9
Accession No.	Uniprot ID: Q9BY32
Uniprot	Q9BY32
GeneID	3704;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ITPA

Note: This product is for in vitro research use only