

## Interleukin-13 receptor subunit alpha-2

Catalog No: #AP84253



Package Size: #AP84253-1 50ug #AP84253-2 100ug #AP84253-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-13 receptor subunit alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	IL-13 receptor subunit alpha-2,IL-13R subunit alpha-2,IL-13R-alpha-2,IL-13RA2,CD213a2
Accession No.	Uniprot ID:Q95LF0
Uniprot	Q95LF0
GeneID	403622;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL13RA2

Note: This product is for in vitro research use only