

Indoleamine 2,3-dioxygenase 1

Catalog No: #AP84254



Package Size: #AP84254-1 50ug #AP84254-2 100ug #AP84254-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Indoleamine 2,3-dioxygenase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IDO-1, Indoleamine-pyrrole 2,3-dioxygenase, Ido, Indo
Accession No.	Uniprot ID:Q9ERD9
Uniprot	Q9ERD9
GeneID	66029;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ido1

Note: This product is for in vitro research use only