

Interleukin-1 receptor-associated kinase 1

Catalog No: #AP84258



Package Size: #AP84258-1 50ug #AP84258-2 100ug #AP84258-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor-associated kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IRAK-1,IRAK
Accession No.	Uniprot ID:P51617
Uniprot	P51617
GeneID	3654;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IRAK1

Note: This product is for in vitro research use only