

ATP-sensitive inward rectifier potassium channel 10

Catalog No: #AP84259



Package Size: #AP84259-1 50ug #AP84259-2 100ug #AP84259-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATP-sensitive inward rectifier potassium channel 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Inward rectifier K(+) channel Kir4.1, Potassium channel, inwardly rectifying subfamily J member 10
Accession No.	Uniprot ID:Q9JM63
Uniprot	Q9JM63
GeneID	16513;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Kcnj10

Note: This product is for in vitro research use only