

Transcription factor AP-1

Catalog No: #AP84266



Package Size: #AP84266-1 50ug #AP84266-2 100ug #AP84266-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Transcription factor AP-1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | AH119,Activator protein 1,AP1,Proto-oncogene c-Jun,V-jun avian sarcoma virus 17 oncogene homolog,Jun A |
| Accession No. | Uniprot ID:P05627 |
| Uniprot | P05627 |
| GeneID | 16476; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Jun

Note: This product is for in vitro research use only