

Hyaluronan synthase 2

Catalog No: #AP84268



Package Size: #AP84268-1 50ug #AP84268-2 100ug #AP84268-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hyaluronan synthase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Hyaluronate synthase 2,Hyaluronic acid synthase 2,HA synthase 2
Accession No.	Uniprot ID:Q92819
Uniprot	Q92819
GeneID	3037;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HAS2

Note: This product is for in vitro research use only