

Hexokinase-3

Catalog No: #AP84269



Package Size: #AP84269-1 50ug #AP84269-2 100ug #AP84269-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hexokinase-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hexokinase type III, HK III
Accession No.	Uniprot ID:Q3TRM8
Uniprot	Q3TRM8
GeneID	212032;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hk3

Note: This product is for in vitro research use only