

Keratin, type I cytoskeletal 19

Catalog No: #AP84270



Package Size: #AP84270-1 50ug #AP84270-2 100ug #AP84270-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 19
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytokeratin-19,CK-19,Keratin-19,K19,Krt1-19
Accession No.	Uniprot ID:P19001
Uniprot	P19001
GeneID	16669;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Krt19

Note: This product is for in vitro research use only