

cAMP-dependent protein kinase catalytic subunit beta

Catalog No: #AP84271



Package Size: #AP84271-1 50ug #AP84271-2 100ug #AP84271-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	cAMP-dependent protein kinase catalytic subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PKA C-beta, Pkacb
Accession No.	Uniprot ID:P68181
Uniprot	P68181
GeneID	18749;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Prkacb

Note: This product is for in vitro research use only