

# Interferon-induced helicase C domain-containing protein 1

Catalog No: #AP84273

Package Size: #AP84273-1 50ug #AP84273-2 100ug #AP84273-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interferon-induced helicase C domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Clinically amyopathic dermatomyositis autoantigen 140 kDa,CADM-140 autoantigen,Helicase with 2 CARD domains,Helicard,Interferon-induced with helicase C domain protein 1,Melanoma differentiation-associated protein 5,MDA-5,Murabutide down-regulated protein,RIG-I-like receptor 2,RLR-2,RNA helicase-DEAD box protein 116,MDA5,RH116
Accession No.	Uniprot ID:Q9BYX4
Uniprot	Q9BYX4
GeneID	64135;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IFIH1

Note: This product is for in vitro research use only