

# Interferon-induced helicase C domain-containing protein 1

Catalog No: #AP84274

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP84274-1 50ug #AP84274-2 100ug #AP84274-3 1mg

## Description

Product Name	Interferon-induced helicase C domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Helicase with 2 CARD domains,Helicard,Interferon induced with helicase C domain protein 1,Melanoma differentiation-associated protein 5,MDA-5,RIG-I-like receptor 2,RLR-2
Accession No.	Uniprot ID:Q8R5F7
Uniprot	Q8R5F7
GeneID	71586;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ifih1

Note: This product is for in vitro research use only