

Interferon-induced transmembrane protein 3

Catalog No: #AP84275



Package Size: #AP84275-1 50ug #AP84275-2 100ug #AP84275-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon-induced transmembrane protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dispanin subfamily A member 2b,DSPA2b,Fragilis protein,Interferon-inducible protein 15,Mouse ifitm-like protein 1,Mil-1
Accession No.	Uniprot ID:Q9CQW9
Uniprot	Q9CQW9
GeneID	66141;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ifitm3

Note: This product is for in vitro research use only