

Glutamine synthetase

Catalog No: #AP84279



Package Size: #AP84279-1 50ug #AP84279-2 100ug #AP84279-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamine synthetase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	GS, Glutamate decarboxylase, Glutamate--ammonia ligase
Accession No.	Uniprot ID:Q8HZM5
Uniprot	Q8HZM5
GeneID	403443;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GLUL

Note: This product is for in vitro research use only