

Keratin, type I cytoskeletal 16

Catalog No: #AP84284



Package Size: #AP84284-1 50ug #AP84284-2 100ug #AP84284-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytokeratin-16,CK-16,Keratin-16,K16,KRT16A
Accession No.	Uniprot ID:P08779
Uniprot	P08779
GeneID	3868;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT16

Note: This product is for in vitro research use only