

Kynurenine--oxoglutarate transaminase 1

Catalog No: #AP84285



Package Size: #AP84285-1 50ug #AP84285-2 100ug #AP84285-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kynurenine--oxoglutarate transaminase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cysteine-S-conjugate beta-lyase,Glutamine transaminase K,GTK,Glutamine--phenylpyruvate transaminase,Kynurenine aminotransferase 1,Kynurenine aminotransferase I,KATI,Kynurenine--oxoglutarate transaminase I,CCBL1
Accession No.	Uniprot ID:Q16773
Uniprot	Q16773
GeneID	883;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KYAT1

Note: This product is for in vitro research use only