

Mast/stem cell growth factor receptor Kit

Catalog No: #AP84286



Package Size: #AP84286-1 50ug #AP84286-2 100ug #AP84286-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mast/stem cell growth factor receptor Kit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	SCFR,Proto-oncogene c-Kit,Tyrosine-protein kinase Kit,CD117
Accession No.	Uniprot ID:Q2HWD6
Uniprot	Q2HWD6
GeneID	396810;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KIT

Note: This product is for in vitro research use only