

Glandular kallikrein-3, submandibular

Catalog No: #AP84287

Package Size: #AP84287-1 50ug #AP84287-2 100ug #AP84287-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glandular kallikrein-3, submandibular
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	rGK-3,RSGK-50,S1 kallikrein,Tissue kallikrein,Klk-3
Accession No.	Uniprot ID:P15950
Uniprot	P15950
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Klk3

Note: This product is for in vitro research use only