

## Keratin, type I cytoskeletal 9

Catalog No: #AP84288

Package Size: #AP84288-1 50ug #AP84288-2 100ug #AP84288-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Keratin, type I cytoskeletal 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokeratin-9,CK-9,Keratin-9,K9,Spermatid perinuclear ring manchette protein K9,Krt1-9
Accession No.	Uniprot ID:Q8CIS9
Uniprot	Q8CIS9
GeneID	266717;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Krt9

Note: This product is for in vitro research use only