

Keratin, type II cytoskeletal 7

Catalog No: #AP84289



Package Size: #AP84289-1 50ug #AP84289-2 100ug #AP84289-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type II cytoskeletal 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytokeratin-7,CK-7,Keratin-7,K7,Sarcolectin,Type-II keratin Kb7,SCL
Accession No.	Uniprot ID:P08729
Uniprot	P08729
GeneID	3855;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT7

Note: This product is for in vitro research use only