

# Eukaryotic translation initiation factor 5A-1

Catalog No: #AP84291



Package Size: #AP84291-1 50ug #AP84291-2 100ug #AP84291-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Eukaryotic translation initiation factor 5A-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	eIF-5A-1,Eukaryotic initiation factor 5A isoform 1,eIF-5A,eIF-4D
Accession No.	Uniprot ID:Q3T1J1
Uniprot	Q3T1J1
GeneID	287444;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Eif5a

Note: This product is for in vitro research use only