

Interferon epsilon

Catalog No: #AP84292



Package Size: #AP84292-1 50ug #AP84292-2 100ug #AP84292-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon epsilon
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IFN-epsilon,Interferon epsilon-1,Interferon tau-1,Ifne1,Ifnt1
Accession No.	Uniprot ID:Q80ZF2
Uniprot	Q80ZF2
GeneID	230405;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ifne

Note: This product is for in vitro research use only